


## EV-300 SEMI RIGID LUMBAR BACK SUPPORT

### Features:

Lumbosacral orthosis fabricated from dense elastic material providing adjustable compression for greater lumbar containment. Incorporating a semi-rigid posterior with four malleable steel stays, threadless seam technology and a removable lumbar cushion. Also featuring a specifically designed fastening system for ease of fitting for patients with limited dexterity.

### Indications:

- Lumbago.
- Lumbosciatica.
- Degenerative processes.
- Weak and atonic muscles.
- Pre- and post-surgical treatment.
- Disc hernia.
- Spondylarthrosis.
- Spondylolysis.



BLACK	BEIGE	SIZE	SIZE	WAIST PERIMETER IN CM
EV-300 XXS	EV-200 XXS	1	XX-SMALL	63 - 75
EV-300 XS	EV-200 XS	2	X-SMALL	75 - 83
EV-300 SM	EV-200 SM	3	SMALL	83 - 91
EV-300 MD	EV-200 MD	4	MEDIUM	91 - 99
EV-300 LG	EV-200 LG	5	LARGE	99 - 107
EV-300 XL	EV-200 XL	6	X-LARGE	107 - 115
EV-300 XXL	EV-200 XXL	7	XX-LARGE	115 - 123

