



EXO-LITE UNIVERSAL BACK BRACE

Features:

The Exo-Lite Universal Back Brace support is a new lightweight low profile option that provides stability and support for chronic sprains and strains of the lower spine, spinal stenosis, and post-surgery support. The mesh material allows for greater breathability of the brace.

It also provides additional support to the lower spine to prevent future injury and will achieve even greater compression for the patient.

- Advanced pulley system allows for tailored compression.
- Mesh material allows for more breathability and a superior low profile, lightweight design.
- Simple strap design allows for a universal fit.

Indications:

- Mild lower back pain.
- Moderate lower back strains and sprains.
- Ideal for pulled or weak muscles.



ITEM	SIZE	WAIST CIRCUMFERENCE INCHES
76-210	UNIVERSAL	25" - 60"

