



UNIVERSAL FOAM KNEE IMMOBILIZER

Features:

This lightweight, foam knee immobilizer provides a wide range of fittings. A key benefit is the reduction of migration due to the five elastic circumferential straps holding it firmly in place. Malleable medial and lateral stays provide customization for additional support. Sturdy posterior stays keep the leg in full extension.



ITEM	SIZE	THIGH CIRCUMFERENCE
223690	12" Universal	Up to 20"
223660	16" Universal	Up to 24½"
223665	18" Universal	Up to 26½"
223670	20" Universal	Up to 29"
223677	20" Universal	Up to 33"
223675	22" Universal	Up to 31½"
223680	24" Universal	Up to 33¼"
223687	24" Universal	Up to 37"