

FP-74 MANUTEC BREATHABLE THUMB IMMOBILIZING SPLINT

Features:

Immobilizing splint in thermoplastic which includes a very comfortable breathable interface. The splint has a small micro-velcro fastener which can be adjusted to fit to the required length.

Function:

Immobilization of the metacarpophalangeal joints of the thumb (post-trauma and post-surgical).

Indications:

- Tendonitis
- Painful and inflammatory episodes
- Instabilities of metacarpophalangeal joints
- Sports injuries (Eg. Skier's thumb)
- Irritations of the metacarpophalangeal joints
- Immobilization of the trapezium-metacarpal joints (Rhizarthrosis of the thumb)
- Immobilization of the metacarpophalangeal joints in the thumb (after injury or surgery)



Replacement Liners are available

LEFT	RIGHT	SIZE	WRIST CIRCUMFERENCE IN CM	LENGTH IN CM
FP-I74-SM	FP-D74-SM	SMALL	<17	11
FP-I74-MD	FP-D74-MD	MEDIUM	17 - 19	12
FP-I74-LG	FP-D74-LG	LARGE	>19	13

