

## FSX AIR WALKER



- Durable, semi-rigid shell helps support the limb while providing full-shell protection
- Integrated inflation system, no need for a separate hand bulb or turning of a valve
- Low profile rocker sole enables a natural gait pattern
- Separate medial and lateral aircell bladder for customized inflation to accommodate changes in swelling patterns
- Covered toe pad for additional protection



**Available in High Top or Low Top versions**

SIZE	HIGH TOP	LOW TOP	MEN'S SHOE SIZE	WOMEN'S SHOE SIZE
SMALL	T98422	T98522	4 - 6	5 - 7
MEDIUM	T98424	T98524	7 - 10	8 - 11
LARGE	T98426	T98526	11 - 13	12 - 14

### Indications:

- Ankle sprains and strains;
- Forefoot sprains;
- Stable foot fractures;
- Distal tibial fractures;
- Distal fibular fractures;
- Syndesmosis (High Ankle Sprains)