



Gripper™ 3 12" Hinged Knee Brace



The Gripper™ 3 Hinged Knee Brace offers you the ability to provide your patients with a custom fitted device that delivers superior medial & lateral stability and resistance to distal migration.



Enhanced Support & Custom Fit

Multiple condylar pads may be repositioned medially or laterally to provide enhanced support and a custom fit for patients with a valgus or varus knee. Easy-peel outer closure straps may also be trimmed to provide a custom fit.

Compliance & Migration

Internal SkinLoc™ pad works in concert with the outer closure strap to grip the proximal end of the calf muscle and inhibit the brace from migrating distally, thereby enhancing patient compliance.

Durability

Metal hinge pockets are reinforced on the inside with a woven nylon fabric at wear points to prevent hinge from wearing through the pocket.

Ease of Application & Comfort

The wraparound anterior closure allows for ease of application and an optimum fit. Four-way stretch fabric in the popliteal fossa minimizes bunching during flexion to improve comfort. Anterior and sides of brace are made with CoolFlex™ material that is breathable and supportive.

Indications:

- Meniscal injuries
- MCL or LCL injuries or sprains
- Medial or lateral knee instabilities

Suggested HCPCS:
L1820 or L1821

Item #	Size	6" Above	Mid-Patella	6" Below
117231	X-Small	11" - 13"	10" - 12"	10" - 12"
117232	Small	13" - 16"	12" - 14"	12" - 14"
117234	Medium	16" - 18"	14" - 16"	13" - 15"
117235	Large	18" - 21"	16" - 18"	15" - 17"
117236	X-Large	21" - 24"	18" - 20"	17" - 19"
117237	XX-Large	24" - 27"	20" - 22"	19" - 21"
117238	3X-Large	27" - 31"	22" - 24"	21" - 23"
117239	4X-Large	31" - 35"	24" - 26"	23" - 26"

*Use mid-patella measurement for primary sizing info

Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com

©2024 Medical Specialties Inc.

