

HV-33 HALLUX-VALGUS IN THERMOPLASTIC



Features:

Night-time corrector for bunions, manufactured in thermoplastic and velour, with micro-fastener seal. It includes a cushioning pad on the big toe which, together with the pressure strap, at the height of the bunion, maintains the first metatarsus aligned and in correct position during night-time use, reducing pain and preventing its possible evolution.



Indications:

- Hallux-Valgus correction



RIGHT	LEFT	SIZE	EUROPEAN FOOT SIZE	WOMEN'S US FOOT SIZE	MEN'S US FOOT SIZE
HV-33D-SM	HV-33I-SM	Small / 1	34 - 36	4 - 6	3 - 4.5
HV-33D-MD	HV-33I-MD	Medium / 2	37 - 40	6.5 - 9	5 - 7.5
HV-33D-LG	HV-33I-LG	Large / 3	41 - 46	9.5 - 13	8 - 12