




K80032 POST-OP ADJUSTABLE KNEE BRACE

Designed to provide controlled range of motion for patients recovering from knee surgery or those who have knee injuries or instabilities.



- Comfortable and extremely lightweight, simple to apply;
- Telescoping bars are easily adjusted for a personalized fit;
- Universal brace so there are less SKUs to stock in inventory;
- Trimmable straps for size customization;
- Drop lock features allows brace to be locked out with the push of a button;
- Easy to use ROM hinge with quick adjusting flexion and extension stops to control and limit knee flexion and extension;
- Quick clip buckles allow for easy application and removal of brace;
- **Product Code K80032 - One Universal size**