

NeuroFlex® Restorative™ Kentucky Kollar

Treats conditions such as Torticollis/wry neck.

Features and Benefits:

- Comfortably and safely correct side leaning of the head or neck to the left or right.
- Design allows tension to relax.
- Remold to gradually correct neck posture.
- Forehead strap available on request.
- Adult and pediatric sizing available.



**Pediatric Model
showing optional
headband**



Adult Model

Kentucky Kollar	Sizing	Forehead to Shoulder
X-Large 20735-XL	Adult	7" - 9"
Large 20735-L	Adult	6" - 7"
Medium 20735-M	Youth	5" - 6"
Small 20735-S	Child	4" - 5"
X-Small 20735-XS	Infant	3.5" - 4.5"
XX-Small 20735-XXS	Newborn	2.5" - 3.5"
Kentucky Kollar Headband	20735 - HB	Universal

The NeuroFlex® Restorative™ Kentucky Kollar is bilateral and fits left or right side.



Phone: 1-866-334-4441
 Fax: 1-888-484-4431
 Email: csr@poportho.com
 Online: www.poportho.com



Professional Orthopedic Products