

MANUTEC FIX IMMOBILIZING WRIST WITH PALM SPLINT



A new closure and adjustment system consisting of a forearm strap and wrist closure, whose single strap design enables wrist and hand adjustment.

In the palmar area, a double Velcro strap system enables adjustment to the width of the hand's metacarpal area.

Characteristics:

Wrist splint manufactured from breathable velour, made using the latest generation materials (triple layer: foam, velour, sponge), with a triple immobilization system comprising a malleable aluminium palm splint and medial-lateral anatomically shaped strips for adaptation.

Indications:

- Wrist immobilization.
- Post-trauma.
- Post-surgery.
- Painful or swollen sequelae.
- Instability.
- Rehabilitation.

ITEM	SIZE	WRIST CIRCUMFERENCE IN CM	WRIST
M760-SM	Small/Medium	14 - 18	
M760-LG	Large/X-Large	18 - 22	

