

## M770 BREATHABLE THUMB IMMOBILIZING SPLINT



### Features:

The splint is made of breathable velour, which has venting holes on the surface to facilitate the movement of air. The immobilizer system consists of a malleable aluminum splint which can be shaped to the morphology of the patient.

The bilateral design can be adapted for both the left and right hand.

The closing system features a velour strap, with hook and loop closure.



### Indications:

- Tendonitis;
- Painful or swollen joints;
- Instability of the first CMC joint;
- Injuries (eg.: Skier's Thumb, Gamekeeper's Thumb);
- Irritations of the first CMC joint;
- Immobilization of the first CMC joint (Thumb arthrosis).



Bilateral design fits both left and right hand

ITEM	SIZE	WRIST CIRCUMFERENCE IN CM (A)	LENGTH IN CM (B)
M770-SM	SMALL	14 - 18	11
M770-LG	LARGE	18 - 22	13

