



M780 MANUTEC FIX RIZART PLUS CMC THUMB BRACE



Features:

Highly resistant and very soft material. The brace is heat-moldable, flexible and comfortable, and perfectly adapts to the shape of the hand.

The design with double injection of materials and the rounded edges guarantee maximum comfort of use while the fastening system, secured by a metal core, enables the degree of stability of the trapeziometacarpal joint to be adjusted as required.

Function:

Stabilizes the base joint of the thumb and places the first metacarpal in its functional position. Places the thumb in an anatomically correct position and enables the other fingers to move properly, ensuring the functionality of the hand.

Relieves pain and reduces the symptoms of loss of strength.

Indications:

- Osteoarthritis of the first CMC joint.
- Post-surgical treatment of the first CMC joint.
- Arthroplasty of the first CMC joint.
- Instability of the first CMC joint.
- Rhizarthrosis.



SIZE	LEFT	RIGHT	PERIMETER CM
SMALL	M780I-SM	M780D-SM	16 - 19
MEDIUM	M780I-MD	M780D-MD	19 - 23
LARGE	M780I-LG	M780D-LG	23 - 26

