

Manutec[®] Fix Rizart Advance

Immobilising abduction thumb splint



► MAIN FEATURES



Breathable, seamless material with holes on the surface for air passage.



An immobilising system comprised of a double plate with a wrap-around design and thumb abduction.



One-hand quick-close for easier fitting.



Manutec® Fix Rizart Advance

- ▶ Thumb immobilisation splint for restricted abduction.
- ▶ Maintains the wrist and finger movement range and completely frees the base of the thumb to relieve bothersome pressure in this area.
- ▶ Equipped with an immobilisation system consisting of a malleable aluminium core to match the patient's morphology.
- ▶ Bilateral in design, it can be adapted to the left or right hand.

Ref.: M790

Immobilising abduction thumb splint

- ▶ **Effects**
 - Thumb support, stabilisation and immobilisation.
 - Abduction restriction.
- ▶ **Recommended for**
 - Rhizarthrosis and osteoarthritis of the trapezium-metacarpal joint.
 - Thumb tendonitis and sprains (skier's thumb).
 - Painful and inflammatory after-effects.
 - Instability in the trapezium-metacarpal and metacarpophalangeal joints.
 - Post-surgical and post-trauma treatment.
- ▶ **Colour**
 - Grey-Orange.

SIZE	PERIMETER /cm/	LENGTH /cm/
1	14-18	11
2	18-22	13



Bilateral design



Distributed in Canada by:
Professional Orthopedic Products
1-866-334-4441
csr@poportho.com

