

NEW DESIGN

ASO[®] MAX[™]



The Choice of ChampionsSM

With the new improvements to the ASO ankle stabilizer, the ASO maintains its leadership position in the ankle brace market. Give your athletes and patients the patented and proven technology of the ASO ankle stabilizer. The leader continues to stay several steps ahead...



Stabilizing Straps

The stabilizing straps are now angled to better contour above the malleoli for enhanced comfort and better migration resistance. Nylon stabilizing straps form a complete figure-eight to help stabilize the subtalar and talocrural joints.

Elastic Cuff

Circumferential woven elastic cuff keeps laces and stabilizing straps secure. In addition, it works in concert with the upper boot to resist the inversion or eversion forces.

CoolFlex Achilles Pad & Tongue

CoolFlex[™] padding has been increased on the inside of the boot along the Achilles tendon to provide greater comfort. The CoolFlex tongue has been contoured to minimize bunching and allow for neater stitching.

Versatility

With the addition of Skinloc[™] lining, the ASO[®] MAX[™] can be worn with or without a sock.

Indications:

- Inversion or Eversion Ankle Sprains
- Medial or Lateral Ankle Instability
- Subtalar Joint Instability

Bilateral design fits left or right foot

Product Code	Size	Ankle Circum.*
264041	X-Small	10" - 11"
264042	Small	11" - 12"
264044	Medium	12" - 13"
264045	Large	13" - 14"
264046	X-Large	14" - 15"
264047	XX-Large	15" - 16"

*measure ankle circumference according to diagram below

Men's Shoe Size

6-7	7 1/2-9	9 1/2-11	11 1/2-13	13 1/2-15	15+
Narrow Foot	pedi	xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	xxl
					xxxl

Women's Shoe Size

5-6	6 1/2-8	8 1/2-10	10 1/2-12	12 1/2-14
Narrow Foot	pedi	xs	s	m
Regular Foot	xs	s	m	l
Wide Foot	s	m	l	xl



Phone: 1-866-334-4441
Fax: 1-888-484-4431



Online: www.poportho.com
Email: csr@poportho.com

