



## MF 93 SEMI-RIGID WRIST SUPPORT WITH PALMAR/DORSAL SPLINTS 2ND AND 3RD METACARPALS

### Features:

The MF-93 features a breathable, lightweight fabric base with a layer of cotton in contact with the skin providing greater comfort and enabling free movement of the fingers.

The palmar and dorsal splints are made from malleable aluminium for effective and separate bracing of the wrist and fingers according to the desired angles.

Fitting is simple thanks to the radial elastic fabric and closure is fast by means of trimmable Velcro straps and anti-torsion buckles.

### Function:

Immobilization of the wrist, hand, 2nd or 3rd fingers. Pain relief.

### Indications:

- Fracture of the 2nd or 3rd metacarpals (boxer's fracture).
- Injuries to the wrist, hand and 2nd or 3rd fingers.
- Dislocations and sprains.
- Pre/post-surgical and post-injury treatment.



LEFT	RIGHT	SIZE	WRIST CIRCUMFERENCE IN CM
MF-I93 SM	MF-D93 SM	SMALL	14 - 16
MF-I93 MD	MF-D93 MD	MEDIUM	16 - 18
MF-I93 LG	MF-D93 LG	LARGE	18 - 20



## MF-95 SEMI-RIGID WRIST SUPPORT WITH PALMAR/DORSAL SPLINTS 4TH AND 5TH METACARPALS

### Features:

The MF-95 features a breathable, lightweight fabric base with a layer of cotton in contact with the skin providing greater comfort and enabling free movement of the fingers.

The palmar and dorsal splints are made from malleable aluminium for effective and separate bracing of the wrist and fingers according to the desired angles.

Fitting is simple thanks to the radial elastic fabric and closure is fast by means of trimmable Velcro straps and anti-torsion buckles.

### Function:

Immobilization of the wrist, hand, 4th or 5th fingers. Pain relief.

### Indications:

- Fracture of the 4th or 5th metacarpals (boxer's fracture).
- Injuries to the wrist, hand and 4th or 5th fingers.
- Dislocations and sprains.
- Pre/post-surgical and post-injury treatment.



LEFT	RIGHT	SIZE	WRIST CIRCUMFERENCE IN CM
MF-I95 SM	MF-D95 SM	SMALL	14 - 16
MF-I95 MD	MF-D95 MD	MEDIUM	16 - 18
MF-I95 LG	MF-D95 LG	LARGE	18 - 20