

## MF-52 SEMI RIGID WRIST SUPPORT WITH PALMER SPLINT



### Features:

The MF-52 Semi Rigid Wrist Support with Palmer Splint features a breathable, lightweight fabric base with a layer of cotton in contact with the skin providing greater comfort and enabling free movement of the fingers. The palmar splint is made from malleable aluminum and it features 2 plastic dorsal supports for effective bracing of the wrist according to the desired angle. Fitting is simple thanks to the dorsal elastic fabric. Fastening the brace can be done quickly by means of the trimmable Velcro straps and anti-torsion buckles.

### Function:


Functional immobilization of the wrist. Stabilization of the carpal joints. Pain relief.

### Indications:

- Carpal tunnel syndrome.
- Wrist sprains.
- Osteoarthritis.
- Tendinitis.
- Overloading of the wrist joints.
- Residual instability.
- Pre/post-surgical and post-injury treatment.



RIGHT	LEFT	SIZE	PERIMETER IN CM	LENGTH IN CM
MF-D52 SM	MF-I52 SM	Small / 1	12 - 16	17
MF-D52 MD	MF-I52 MD	Medium / 2	16 - 18	17
MF-D52 LG	MF-I52 LG	Large / 3	18 - 20	17



## MP-72 SEMI RIGID WRIST SUPPORT WITH THUMB SPLINT



### Features:

The MP-72 Semi Rigid Wrist Support with Thumb Splint features a breathable, lightweight fabric base with a layer of cotton in contact with the skin providing greater comfort and enabling free movement of the fingers. The thumb splint is made from malleable aluminum and it features 2 plastic dorsal supports for effective bracing of the thumb according to the desired angles. In addition, the radial splint features a semi-circular design for increased comfort and improved positioning of the thumb. Fitting is simple thanks to the dorsal elastic fabric and closure is fast by means of trimmable Velcro straps and anti-torsion buckles.



### Function:

Functional immobilization of the thumb. Stabilization of the carpal joints. Pain relief.

### Indications:

- Rhizarthrosis.
- Ulnar collateral ligament injury of the thumb (skier's thumb).
- Tendinitis.
- De Quervain's tenosynovitis
- Overloading of the wrist and thumb joints.
- Residual instability.
- Pre/post-surgical and post-injury treatment.



RIGHT	LEFT	SIZE	PERIMETER IN CM	LENGTH IN CM
MP-D72 SM	MP-I72 SM	Small / 1	12 - 16	21
MP-D72 MD	MP-I72 MD	Medium / 2	16 - 18	21
MP-D72 LG	MP-I72 LG	Large / 3	18 - 20	21

