

# Wrist Supports



## Motion Manager<sup>®</sup> Wrist Support

This neoprene wrist support is an excellent choice for providing functional motion. The contoured plastic palmar stay and the multi-density dorsal stays offer adjustability to provide the optimum level of flexion and extension control.

A low profile design allows the range of motion needed to return to work safely.

Bilateral design fits left or right hand.



Item #	Size	Wrist Circ.
223901	X-Small	up to 5 <sup>3</sup> / <sub>4</sub> "
223902	Small	5 <sup>3</sup> / <sub>4</sub> " - 6 <sup>1</sup> / <sub>2</sub> "
223904	Medium	6 <sup>1</sup> / <sub>2</sub> " - 7 <sup>1</sup> / <sub>4</sub> "
223905	Large	7 <sup>1</sup> / <sub>4</sub> " - 8"
223906	X-Large	over 8"

Retail Packaging

U.S. Patent #5,928,172

Phone: 1-866-334-4441  
Fax: 1-888-484-4431



Online: [www.poportho.com](http://www.poportho.com)  
Email: [csr@poportho.com](mailto:csr@poportho.com)