



NAVIGAIT

Features:

NAVIGAIT™ is an alternative for users with mild foot drop. The amount of dorsiflexion is easily adjusted with an elastic coil to accommodate the patient's specific needs.

Supports ankle dorsiflexion during swing phase. Elastic coil adjustable and easily replaceable. Contains removable interface. The amount of dorsiflexion is adjusted through the position in the lacing and with the velcro® attachment on the shin plate. Low profile design makes it invisible under slacks. Anterior shell can be trimmed with scissors if needed.

The optional 4-FOOT Wrap accessory makes it possible to wear NAVIGAIT™ with sandals, shoes that do not have an anchor point, or barefoot. The wrap has a plastic insert for increased stabilization and non-slip silicon strips to keep it in place.



Recommended Range Of Application:

Intended to support a foot when the ability to actively dorsiflex is reduced due to conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, or Spinal Cord Injury.

Contraindications:

NAVIGAIT™ should not be used by patients who present with spasticity or severe edema.

NAVIGAIT

ITEM	SIZE	LENGTH BETWEEN KNEE AND ANKLE
284100010	XS	UNDER 9"
284100011	SM/MD	UNDER 11"
284100013	LG/XL	OVER 11"

4-FOOT WRAP

ITEM	SIZE	FORE-FOOT CIRCUMFERENCE
284150011	SMALL	8.25" - 9"
284150012	MEDIUM	9" - 10.25"
284150013	LARGE	10.25" - 11.25"

