



NEUROFLEX™ RESTORATIVE™ FLEX BOOT

Function:

Safely and comfortably treats inversion or eversion, foot/ankle contractures, foot drop and plantar flexion.

Features and Benefits:

- Suspends the heel.
- Heat mold to treat almost any foot/ankle deformity.
- Remold to gradually bring joint toward and to correction.
- Semi-rigid design with flex properties allows the joints to reach relaxation then continues the slow gentle stretch.
- Non-Skid sole upon request.



PRODUCT	SIZE	FOOT LENGTH
20734-S	SMALL	Up to 8"
20734-M	MEDIUM	Up to 9½"
20734-L	LARGE	Up to 11"