



## NEUROFLEX™ RESTORATIVE™ FLEX KNEE

### Function:

Provides prolonged low load passive stretch and treats flexion contractures of the knee that result in lost range of motion. Addresses mild to severe flexion contractures.

### Features and Benefits:

- Padding provides comfort and reduces skin pressure.
- High-temperature plastic cuffs can be remolded by your health professional to custom fit the patient's upper and lower legs.
- Adjustable dials allow gradual changes to move the joint toward normal alignment. Locks in to apply mild stretch. Gradually reset dials to work joint toward normal alignment and re-lengthen shortened tissue.
- Brace provides 3-point leverage similar to manual stretching, improving range of motion.
- Flex Technology splint moves with the patient's abnormal muscle tone/spasticity. Side bars flex when patient draws inward, then brings the joint back to preset position.
- Comfortable.



PRODUCT	SIZE	TTK	KTC	THIGH CIRC.	CALF CIRC.
20725C	CHILD	5"	6"	8" - 12"	7" - 10"
20725S	SMALL	6½"	6½"	10" - 15"	8" - 12"
20725M	MEDIUM	7½"	7½"	13" - 17"	10" - 15"
20725L	LARGE	8½"	8½"	16" - 22"	14" - 18"

