


Neuroflex™ Restorative™ Elbow

Treats mild to severe contractures of the arm.

Features:

- Restore normal range of motion with prolonged low load passive stretch, with flex properties to relax joint.
- Adjustable dials lock in to apply mild stretch at elbow joint to achieve a relaxation.
- Side bars flex when patient draws inward, then brings the joint back to preset position.
- Gradually reset dials to work joint toward normal alignment and re-lengthen shortened tissue.
- Comfortable.



ELBOW MEASURING GUIDE

Size	BTE	ETF	Biceps Cir.	Forearm Cir.
Child	3½"	4"	5½-7½"	3½-6½"
Youth	4"	5"	6½-8½"	5½-7½"
Ex-Small*	5"	5"	7½-9½"	6½-8½"
Small	5"	6"	8-12"	7-10"
Medium	6"	6"	8-13"	7½-11"
Large	6½"	6½"	10-15"	8-12"

*Ex-small is a new size to help bridge the pediatric and adult sizing

Phone: 1-866-334-4441
Fax: 1-888-484-4431



Online: www.poportho.com
Email: csr@poportho.com