

O-104 NAVEL TRUSS


Features:

Elastic strap with rear Velcro fastening and interchangeable navel cushion.

Indications:

- Reduction of umbilical hernias.
- Good compression on the navel area.

ITEM	SIZE	WAIST CIRCUMFERENCE IN CM
O-104-7	7	73 - 83
O-104-8	8	84 - 91
O-104-9	9	92 - 99
O-104-10	10	100 - 107
O-104-11	11	108 - 120



O-105 NAVEL TRUSS

Features:

Elastic belt made of breathable, compressive fabric with adjustable closure. It has a movable silicone pad lined with cotton fabric. With semi-rigid reinforcements for proper adaptation and to prevent it from folding.

Indications:

- Reduction of umbilical hernias.
- Good compression on the navel area.

ITEM	SIZE	WAIST CIRCUMFERENCE IN CM
O-105-1	1	70 - 95
O-105-2	2	95 - 125
O-105-3	3	125 - 150

