



ELASTIC KNEE SUPPORT WITH LATERAL STABILIZERS



Features:

Manufactured in a breathable stretch knit fabric that is very strong yet so soft for added comfort; its seamless construction means no uncomfortable rubbing. Three dimensional anatomical adjustment. The product includes a ring-shaped viscoelastic pad and lateral stabilizers for improved medial-lateral control.

Function:

The massaging effect of the compressive fabric and knee pad improves circulation and promotes faster reabsorption of bruises and swelling. Comfort zones made of soft knit fabric avoid pressure and undesired compression.



Indications:

- Osteoarthritis
- Arthrosis
- Minor knee injuries and instability
- Inflammatory processes
- Support and protection for athletic and work activities

Bilateral design fits both left and right knee

24 CM LENGTH	28 CM LENGTH	SIZE	KNEE CIRCUMFERENCE IN CM
OS6211-SM	OS6212-SM	SMALL	32 - 36
OS6211-MD	OS6212-MD	MEDIUM	36 - 41
OS6211-LG	OS6212-LG	LARGE	41 - 48

