



OS6230 ELASTIC ELBOW SUPPORT WITH GEL PADS



Features:

Manufactured in a breathable stretch knit fabric that is very strong yet so soft for added comfort; its seamless construction means no uncomfortable rubbing. Three dimensional anatomical adjustment. The product includes viscoelastic pads that redistribute the compression on the medial and lateral epicondyles, and an elastic strap for individualized adjustment.

Function:

The controlled compression provides pain relief and helps eliminate tension. The pads redistribute the compression from the prominent epicondyles to the surrounding soft tissues and provide an intermittent massaging effect. Around the elbow, the knit fabric is softer, avoiding pressure and undesired compression in this area.

Indications:

- Lateral Epicondylitis (tennis elbow)
- Medial Epicondylitis (golfer's elbow)
- Tendonitis
- Mild contusions
- Inflammatory processes
- Injuries and instability
- Support and protection for athletic and work activities



Bilateral design fits both left and right elbow

ITEM	SIZE	ELBOW CIRCUMFERENCE IN CM
OS6230-SM	SMALL	22 - 26
OS6230-MD	MEDIUM	26 - 30
OS6230-LG	LARGE	30 - 34

