



## OS6240 ELASTIC ANKLE SUPPORT WITH GEL PADS



### Features:

Manufactured in a breathable stretch knit fabric that is very strong yet so soft for added comfort; its seamless construction means no uncomfortable rubbing. Three dimensional anatomical adjustment. The product includes viscoelastic pads that redistribute compression from the prominent malleolus to the surrounding soft tissues.

### Function:

The massaging effect of the compressive fabric and padding stimulates faster reabsorption of bruises and swelling. Around the heel, the knit fabric is softer, avoiding pressure and undesired compression in this area.

### Indications:

- Mild ankle instability
- Joint effusion and inflammation of traumatic origin as well as that caused by previous diseases
- Postoperative and post-injury (after sprains)
- Tendomyopathies
- Support and protection for athletic and work activities



Bilateral design fits both left and right ankle

ITEM	SIZE	ANKLE CIRCUMFERENCE IN CM
OS6240-SM	SMALL	17 - 20
OS6240-MD	MEDIUM	20 - 23
OS6240-LG	LARGE	23 - 26

