



## Pavlik Harness for treatment of Developmental Hip Dysplasia

The Deluxe Pavlik Harness is used in the treatment of Developmental Hip Dysplasia (DHD). DHD is treated by correctly positioning the head of the femur (thigh bone) within the hip socket. This will allow the hip joint to develop normally. The correct position is achieved when the baby's legs are flexed to approximately 90 degrees to 100 degrees and the legs are brought out to the side (abduction).

### Features:

- Colour-coded straps for easy placement;
- Soft plush materials;
- Well padded;
- Machine washable Plastic buckles;
- Durable hook and loop closures;
- Easily changed from closed to open toe;
- Holds the legs in “Frog” position while allowing hip movement ;
- Latex Free.



**INTEGRATED SOCKS**

ITEM	SIZE	CHEST CIRCUMFERENCE IN INCHES
PAVLIK XS	X-SMALL	12” - 14½”
PAVLIK SM	SMALL	14” - 16½”
PAVLIK MD	MEDIUM	16½” - 18½”
PAVLIK LG	LARGE	18½” - 21½”

