



POP PADDED ARM SLING

Designed to provide an optimum level of comfort and shoulder stabilization.

Features:

- Soft, padded shoulder strap.
- Thumb loop.
- Stretch-resistant, lightweight, breathable material.
- Latex-free

Benefits:

- Easily adjustable.
- Increase functionality.
- Increased durability.
- High compliance.



ITEM	SIZE	MEASUREMENTS
5552309	UNIVERSAL	UP TO 16"
5552301	XX-SMALL	UP TO 12"
5552302	X-SMALL	UP TO 13"
5552303	SMALL	UP TO 14"
5552304	MEDIUM	UP TO 15"
5552305	LARGE	UP TO 16"
5552306	X-LARGE	UP TO 17"

To determine proper size, measure the linear distance from the olecranon process (tip of elbow) to the 5th MP joint (knuckles).

