



BENDEASE HAND

Provides moderate, comfortable support of the hand and wrist, post injury or surgery.

Features:

- Foam over aluminum base can be remolded by health professional to continue restorative process as the patient improves.
- Lightweight support for painful extremities.
- Base extends under the thumb gradually abducting it comfortably.
- Padding provides comfort and redistributes skin pressure.
- Color coding aids proper strap attachment.
- Comfortable thumb support for mild lost range.
- Includes optional finger separators which may be removed when not in use.



LEFT	RIGHT	SIZE	WIDTH OF MCP	WRIST TO FINGERS
20715SL	20715SR	SMALL	UP TO 3"	UP TO 6½"
20715ML	20715MR	MEDIUM	UP TO 3½"	UP TO 7"
20715LL	20715LR	LARGE	UP TO 3¾"	UP TO 8"
20715L	20715R	UNIVERSAL		

