



RESTORATIVE HAND

Provides moderate, comfortable support of the hand and wrist, post injury or surgery.

Features:

- High temperature plastic base can be remolded by health professional to continue restorative process as the patient improves.
- Lightweight support for painful extremities.
- Plastic base extends under the thumb gradually abducting it comfortably.
- Flex Technology splint moves with the patient's abnormal muscle tone and spasticity for comfort, and helps to relax the tone.
- Padding provides comfort and reduces skin pressure.
- Color coding aids proper strap attachment.
- Comfortable thumb support for mild lost range.
- Available with a dorsal attachment in all sizes which is ideal for constrained induced movement therapy.

Indications:

Traumatic brain injury / Stroke / Spinal cord injury / Cerebral Palsy / Multiple Sclerosis / Parkinson's Disease


Product Codes:

Adult Regular: 20736

Pediatric Regular: 20717

Adult with Dorsal: 20736-MOD

Pediatric with Dorsal: 20717-MOD



Adult Sizing

SIZE	WIDTH OF MP	WRIST TO FINGERS
X-SMALL	2¾" - 3"	UP TO 5"
SMALL	UP TO 3"	UP TO 6½"
MEDIUM	UP TO 3½"	UP TO 7"
LARGE	UP TO 3¾"	UP TO 8"
X-LARGE	UP TO 4"	UP TO 9"

Pediatric Sizing

SIZE	WIDTH OF MPC	WRIST TO FINGERS
INFANT	1¼" - 2"	UP TO 2½"
TODDLER JR	1¼" - 2"	UP TO 2½"
TODDLER	2" - 2¼"	UP TO 3"
CHILD	2¼" - 2½"	UP TO 4"
YOUTH	2¾" - 3"	UP TO 5"

