

# Restorative Safe Boot™ II

Comfortably treats heel ulcers, plantar flexion, ankle alignment, and foot drop.

## Features and Benefits:

- The flanged heel on the boot elevates the heel and relieves pressure thereby preventing the development of heel ulcers and allowing for healing of pre-existing ulcers.
- Keeps patient's foot in dorsiflexion thereby preventing drop foot/heel cord contractures that may develop due to prolonged bed rest.
- Rotator bar (kickstand), when extended medially or laterally, inhibits internal or external rotation at the hip and therefore prevents hip contractures.
- Heat moldable frame.
- External cover can be removed for laundering.
- Low profile appearance for easy compliance.



Safe Boot™ II  
with foot extension



Moldable Safe Boot™ II  
Frame

SafeBoot™ II	Length of Foot	Pediatric Safeboot™	Length of Foot
Small 20744-II-S	Up to 8"	20746-3*	3"
Medium 20744-II-M	Up to 9.5"	20746-3.5*	3.5"
Large 20744-II-L	Up to 11"	20746-4*	4"
X-Large 20744-II-XL	Up to 13"	20746-5*	5"
		20746-6	6"
		20746-7	7"
		20746-8	8"

The Safe Boot™ II and Pediatric Safeboot™ are bilateral and fit left or right foot.

\*Flex straps not available on 3", 3.5", 4", and 5" models.

Phone: 1-866-334-4441  
 Fax: 1-888-484-4431  
 Email: [csr@poportho.com](mailto:csr@poportho.com)  
 Online: [www.poportho.com](http://www.poportho.com)



Professional Orthopedic Products