

Thumb Supports



TeePee™ Thumb Protector

This unique configuration with its three aluminum stays forms the “TeePee” that comfortably supports the MCP (metacarpal-phalangeal) joint. Each stay can be shaped and contoured with ease to provide a custom fit. The TeePee does not restrict wrist motion and fits the left or right thumb (bilateral).



PDAC Verified HCPCS: L3923

Sizing for the TeePee
To Determine proper size, measure the linear distance from the dorsal crease of the thumb IP joint to the peak of the radial styloid process.



Size	Item Number	Measurements
x-small	223081	under 3 1/2”
small	223082	3 1/2” - 4 1/8”
medium	223084	4 1/8” - 4 5/8”
large	223085	4 5/8” - 5 1/8”
x-large	223086	over 5 1/8”

CMC Thumb Support

The CMC Thumb Support aids in the treatment of arthritis and instability of the CMC (carpometacarpal) joint. The unique design of the internal padded stay in combination with the outer elastic strap provides focal compression at the base of the thumb and helps stabilize the CMC joint. The CMC Thumb Support provides this functional stabilization without restricting the use of the thumb in pinching or grasping objects.



Suggested HCPCS: L3923

Size	Left	Right	Measurements
x-small	224421	224431	under 5 1/2”
small	224422	224432	5 1/2” - 6 1/4”
medium	224424	224434	6 1/4” - 7”
large	224425	224435	7” - 7 3/4”
x-large	224426	224436	over 7 3/4”

*measure circumference of wrist at narrowest point

Medical Specialties, Inc.
4600-K Lebanon Rd.
Charlotte, NC 28227
p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com
www.medspec.com

©2007-2023 Medical Specialties, Inc.



rev 5/23