



## TELESCOPING RANGE OF MOTION ELBOW BRACE

### Features:

The Telescoping Range of Motion Elbow Brace is a lightweight, malleable to the arm and comfortable-to-wear product that allows the patient to resume normal activities of daily living with complete confidence. In the healing phase, comfort means compliance.

### Clinical Product Advantages:

- Universal in size and comes in Left or Right.
- Easy ROM settings from 0-120 degrees.
- One-touch button enables easy setting of telescoping arms.
- Lightweight aluminum upright.
- Comfortable strap and cuff padding reduces migration.
- Durable D-ring design to help restrict movement.
- Adjustable and trimmable medical grade Velcro straps for optimal sizing.
- Malleable Bicep and forearm cuffs.
- Lateral Condyle padding.
- Optional shoulder strap accessory for patients requiring a “sling” position.

### Indications:

- Range of motion control post injury.
- Stable fractures of the elbow, distal humerus, proximal radius or ulna.
- Tendon and ligament injuries or repairs.
- Collateral ligament reconstructions.
- Elbow dislocations.
- Elbow hyperextension.
- Tennis elbow release.
- Triceps tendon repair.
- Ulna nerve injuries or transposition.



| ITEM    | SIZE            |
|---------|-----------------|
| 18-9000 | UNIVERSAL RIGHT |
| 18-9001 | UNIVERSAL LEFT  |