



Professional Orthopedic Products
Providing Innovative Products
Of Unsurpassed Quality



Tel. (866) 334-4441 • Fax (888) 484-4431
csr@poportho.com • www.poportho.com

Celebrating Over 20 Years
Serving Canadians



TGO490 Functional Elastic Ankle Support

Effects: Ankle stabilization, securing and compression / better proprioception.

Indications: Sprains and contusions / Fibromyalgia / Effusion and swelling caused by arthrosis and arthritis / Post-surgical and post-trauma irritation.

ITEM	SIZE	ANKLE CIRC. CM
TGO490-XS	X-Small / 1	19 – 21
TGO490-SM	Small / 2	21 – 23
TGO490-MD	Medium / 3	23 – 25
TGO490-LG	Large / 4	25 – 27
TGO490-XL	X-Large / 5	27 - 29



TGO484 Functional Elastic Calf Support

Effects: Support and compression for the calf muscles / Lower swelling and oedema in the injured muscles / Better proprioception.

Indications: Muscle strain and contusions / Torn or ruptured muscles / Post-trauma irritation.

ITEM	SIZE	CALF CIRC. CM
TGO484-XS	X-Small / 1	27 – 29
TGO484-SM	Small / 2	29 – 32
TGO484-MD	Medium / 3	32 – 36
TGO484-LG	Large / 4	36 – 40
TGO484-XL	X-Large / 5	40 – 45



TGO480 Functional Elastic Knee Support with Strips

Effects: Knee stabilization, securing and compression / Better proprioception.

Indications: Patellar instability / Effusion and swelling caused by arthrosis and arthritis / Fibromyalgia / Sprains and contusions / Post-surgical and post-trauma irritation.

ITEM	SIZE	THIGH CIRC. CM	CALF CIRC. CM
TGO480-XS	X-Small / 1	41 – 44	31 – 34
TGO480-SM	Small / 2	44 – 47	34 – 37
TGO480-MD	Medium / 3	47 – 50	37 - 40
TGO480-LG	Large / 4	50 – 53	40 - 43
TGO480-XL	X-Large / 5	53 - 56	43 - 46



TGO487 Functional Elastic Knee Support with Joints

Effects: Knee stabilization, securing and compression / Better proprioception.

Indications: Patellar and mediolateral instability / Meniscus injuries / Effusion and swelling caused by arthrosis and arthritis / Fibromyalgia / Sprains and contusions / Post-surgical and post-trauma irritation.

ITEM	SIZE	THIGH CIRC. CM	CALF CIRC. CM
TGO487-XS	X-Small / 1	39 - 42	29 – 32
TGO487-SM	Small / 2	42 – 45	32 – 35
TGO487-MD	Medium / 3	45 – 48	35 – 38
TGO487-LG	Large / 4	48 – 51	38 – 41
TGO487-XL	X-Large / 5	51 - 54	41 – 44



TGO481 Functional Elastic Thigh Support

Effects: Support and compression for the thigh muscles / Lower swelling and oedema in the injured muscles / Better proprioception.

Indications: Muscle strain and contusions / Torn or ruptured muscles / Post-trauma irritation.

ITEM	SIZE	THIGH CIRC. CM
TGO481-XS	X-Small / 1	40 – 44
TGO481-SM	Small / 2	44 – 48
TGO481-MD	Medium / 3	48 – 52
TGO481-LG	Large / 4	52 – 57
TGO481-XL	X-Large / 5	57 - 62



TGO330 Functional Elastic Shoulder Support

Effects: Shoulder stabilization, securing and compression / An anti-gravity effect that relieves pressure on the joints by dissipating muscle tension / Better proprioception.

Indications: Muscle strain and contusions / Torn or ruptured muscles / Sprains and subluxation / Fibromyalgia / Arthrosis / Post-surgical and post-trauma irritation.

ITEM	SIZE	BICEP CIRC. CM
TGO330-XS	X-Small / 1	22 – 26
TGO330-SM	Small / 2	26 – 31
TGO330-MD	Medium / 3	31 – 37
TGO330-LG	Large / 4	37 – 44
TGO330-XL	X-Large / 5	44 - 51



TGO340 Functional Elastic Elbow Support

Effects: Pressure on the elbow joint / Better proprioception.

Indications: Epicondylitis and epitrochleitis / Sprains and contusions / Fibromyalgia / Effusion and swelling caused by arthrosis and arthritis / Post-surgical and post-trauma irritation.

ITEM	SIZE	ELBOW CIRC CM
TGO340-XS	X-Small / 1	21 – 23
TGO340-SM	Small / 2	23 – 25
TGO340-MD	Medium / 3	25 – 27
TGO340-LG	Large / 4	27 – 29
TGO340-XL	X-Large / 5	29 - 31



TGO350 Functional Elastic Wrist Support

Effects: Wrist stabilization, securing and compression / Better proprioception.

Indications: Carpal tunnel syndrome / Effusion and swelling caused by arthrosis and arthritis / Fibromyalgia / Sprains and contusions / Post-surgical and post-trauma irritation.

ITEM	SIZE	WRIST CIRC. CM
TGO350-XS	X-Small / 1	14 – 15
TGO350-SM	Small / 2	15 – 16
TGO350-MD	Medium / 3	16 – 17
TGO350-LG	Large / 4	17 – 19
TGO350-XL	X-Large / 5	19 - 21

Basic Concepts of Thera Go®

FUNCTIONAL Therapy

Finding solutions for new trends in trauma and chronic injury rehabilitation is essential at Orliman®. Several studies have found rehabilitation is much faster with early mobilization methods in addition to achieving a higher degree of patient satisfaction. Functional orthotics are designed to improve neuromuscular function, increase static and dynamic balance, and reduce the possibility of repeat injury.

EFFECTIVE Compression

The multi-directional compression provided prevents liquid build-up, reduces stress on the soft tissues and stimulates proprioceptors, thereby enhancing articular function and muscle co-ordination. On the other hand, it optimises blood circulation and facilitates body temperature regulation.

AESTHETIC Design

The modern design for this range is inspired by the anatomic structures of each body segment. The use of organic shapes provides patients with a feeling of protection and movement.

DYNAMIC Rehabilitation

Prolonged or unnecessary immobilization leads to greater functional restriction and, thus, more rigidity in the injured joint or tissue. Thera Go® was developed to make functional recovery easier for patients. Besides planning in collaboration with healthcare professionals, the use of Thera Go® orthotics ideally prepares patients for their daily life and sports.

COMFORTABLE and RESISTANT Fabric

The fabric is a combination of elastic and non-elastic threads for functional and effective rehabilitation.

- Reinforced ends for anatomic adaptation, thereby expanding the useful life of the products.
- Silicone tips on the inside to prevent the support from moving during activity.
- The seams have been minimized on the inside for comfortable use without bothersome chafing.