



▶ FAST-FIX® ANTI-EQUINUS
QR Video

NEW GENERATION

ANKLE FOOT ORTHOSIS Fast-fix®

Graphene-enhanced
engineering and design

A product focused on patients who need walking assistance **to regain their well-being and quality of life.**

NEW
DESIGN



Tested design



Recyclable



100%
made in Spain



Ankle Foot Orthosis (Rancho los amigos - Leaf Spring AFO)

TP-2102D (Right) disassembled

TP-2102I (Left) disassembled

▪ **Indications:**

- In general, patients affected by a limitation in their dorsiflexion of the foot or equinus foot deformity.
- Hemiplegia due to cerebrovascular accidents (stroke), and impairments associated with flaccid paralysis.
- Retraction of the Achilles tendon.
- Lesions to the external popliteal sciatic nerve.
- This orthosis is particularly contraindicated for patients with deviations in the hindfoot varus or valgus and spastic paralysis.

▪ **Effects:**

- Prevents equinus foot deformity / foot drop.

▪ **Colour:**

Black.

SIZE	FOOTWEAR SIZE	HEIGHT (cm)
1	35-37	29
2	38-40	32
3	41-43	34
4	44-46	37



MAIN FEATURES

