



HUMERUS BRACE IN THERMOPLASTIC

Characteristics:

Thermoshaped, two-piece humerus brace in low-density polyethylene, velour adjustment straps and microhook fastening system. Can be cut according to the desired length, coupled with 1 mm polyethylene tabs. Internal plastazote padding.

Indications:

- Diaphyseal fractures of the humerus.
- Post-surgical treatment of fractures and control of the instability of the fractured area during the healing process.



ITEM	SIZE	PROXIMAL IN CM "A"	DISTAL IN CM "B"	VARIANCE
TP-6400XS	X-SMALL	25	22	±3
TP-6400SM	SMALL	27	24	±3
TP-6400MD	MEDIUM	30	27	±4
TP-6400LG	LARGE	32	29	±4
TP-6400XL	X-LARGE	34	31	±4

