



Tripod Lacer IITM Wrist Support



The Tripod Lacer II wrist support offers enhanced immobilization and ease of application. In addition to flexion-extension control, the Tripod Lacer II resists radial-ulnar deviation and its patented design provides equal tension across all laces.

Immobilization

Ulnar, palmar and dorsal aluminum stays work in concert to restrict radial-ulnar deviation and minimize flexion-extension of the wrist. Thumb strap has additional hook/loop interface for added security. Ulnar and palmar stays may be removed if desired for a progressive treatment protocol.

Lacer II Closure System

Patented lace closure design provides equal tension across all laces to easily generate uniform compression for unparalleled support. The tension is easily adjusted for patients with disproportionate forearm and wrist circumferences.

Optimum Comfort

The polymer felt liner and perforated suede allow moisture to evaporate away from the skin and enhance patient compliance.



Indications:

- Post fracture care
- Post operative care
- Carpal tunnel syndrome
- Wrist sprains or strains

Suggested HCPCS: L3908

Size	8" Left	8" Right	Wrist Circumference
X-Small	224001	224011	4 ³ / ₄ " - 5 ¹ / ₂ "
Small	224002	224012	5 ¹ / ₂ " - 6 ¹ / ₄ "
Medium	224004	224014	6 ¹ / ₄ " - 7"
Large	224005	224015	7" - 7 ³ / ₄ "
X-Large	224006	224016	7 ³ / ₄ " - 8 ¹ / ₂ "

