



Spina® II - Spinal Orthosis

SPINA® II Lightweight Plastic Spinal Orthosis

Traditional plastic chairback orthoses lined with plush soft-foam. The 2" straps are made of double-sided Velcro® webbing, adjustable at both ends which provide a wide range of flexibility.

Men's Model I 200, Women's Model I 201



SPINA® II Knight-Taylor Spinal Orthosis

Traditional plastic Knight-Taylor spinal orthoses lined with a 3/8" loop compatible plush soft-foam. The apron front is the same 3/8" plush-foam material with vertical stays and durable 2" plastic loops. The 2" straps are made of double-sided loop webbing, adjustable at both ends which provide a wide range of flexibility. Adjustable over-the-shoulder straps with foam axillary pads.

Men's Model I 275, Women's Model I 276

Sizing

Sizing below applies to both Spina® II models above. Measurements are for male and female/low profile models.

| SIZE | HIP CIRCUMFERENCE IN INCHES |
|---------|-----------------------------|
| SMALL | 30" - 35" |
| MEDIUM | 36" - 41" |
| LARGE | 42" - 47" |
| X-LARGE | 48" - 52" |



Male Front



Female Front