



## VECTOR ABDUCTED SHOULDER SLING


This sling offers superior comfort and shoulder immobilization with its Vector™ strapping system that eliminates pressure on the neck. The elbow region of the sling envelope is open to minimize pressure on the ulnar nerve. The Vector™ shoulder immobilizer is available for sale in either the left or right shoulder configuration. However, each Vector™ is bilateral and may be modified to fit the left or right shoulder. Pillow is designed to position the arm in approximately 15° of abduction.



### Indications:

- Rotator cuff repair
- Bankhart repair
- Shoulder instability
- Glenohumeral dislocation or subluxation
- Soft tissue shoulder strains or sprains

SIZE	LEFT	RIGHT	MEASUREMENT
Small	183702	183712	Up to 13 ½"
Medium	183704	183714	Up to 14 ½"
Large	183705	183715	Up to 15 ½"



To determine the proper size, measure the linear distance from the olecranon process (tip of the elbow) to the 5th MP joint (knuckles).