

Wrist Lacer™ II

Universal wrist support



With its unique closure design, the Wrist Lacer II provides equal tension across all laces offering supreme comfort, immobilization and ease of application.



10½" wrist & forearm support

Universal Design

Dorsal stay casing is easily repositioned to fit a wide range of patient sizes. In addition, dorsal and palmar stays are removable if desired.

Lacer II Closure System

Unique lace closure design provides equal tension across all laces to easily generate uniform compression for unparalleled support. The tension can be easily adjusted for patients with disproportionate forearm and wrist circumferences.

Optimum Comfort

New and improved polypropylene felt liner now has a softer feel and less stretch for enhanced support. Pad inside stockinette protects tendons by cushioning distal lace at base of thumb.

Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Post cast removal
- Postoperative care



8" wrist support



8" wrist support

Suggested HCPCS: L3908

| Size | 10.5" Left | 10.5" Right | Wrist Circ. |
|-----------|------------|-------------|-------------|
| Universal | 223380 | 223381 | 5 ½" - 8 ½" |

| Size | 8" Left | 8" Right | Wrist Circ. |
|-----------|---------|----------|-------------|
| Universal | 223350 | 223351 | 5 ½" - 8 ½" |

U.S. Patent No. 10,772,753

